Seronorm™ CRP Liquid L-1

SERO

REF

213005

LOT

2103819

2...8 °C

2024-05

12 x 1 ml

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Analytical data assignment from invicon

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EN - invicon

Intended use

To use this control material in compliance with a complete Internal QC according the latest RiLiBÄK law.

Assignment and reassignment of target values

The mean analytical values and acceptable ranges indicated and marked in the table below are assigned or reassigned from replicate analyses obtained through collaboration with a number of independent laboratories under the responsibility of invicon. However variations in the analytical value and acceptable range for each component may vary over time due to differences in reagent and/or calibrator lots as well as the instrument and/or method used. If applicable the RiLiBÄK range (Richtlinie der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2019)) is provided.

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EN - SERO (OV)

Product description

Seronorm™ CRP Liquid is a liquid human-based control serum produced from blood collected from voluntary blood donors. Each unit is separately checked using CE approved tests to be negativ for HBs antigen, HIV p24-antigen and HIV I, II and HCV antibodies. However, since no method can completely exclude the presence of infectious agents, this material should be handled as an ordinary patient sample.

Intended use

Seronorm™ CRP Liquid is intended for use as an assayed quality control serum to monitor the precision and trueness of laboratory measurement procedures.

Storage and stability

This product will be stable until the expiry date when stored unopened at 2-8 °C.

Once opened, CRP is stable for 31 days when stored tightly capped in the original vial at 2-8 °C and for 7 days at 15-25 °C.

Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

Instructions for use

1. Swirl the contents gently to ensure homogeneity and the material has achieved room temperature.
2. Remove the amount of material required for testing.
3. Promptly replace the stopper in the vial and store the remaining control material at 2-8 °C.
4. Any material removed from the vial but not used should be discarded. DO NOT RETURN control material to the original vial once removed.

Assignment of values

Seronorm™ CRP Liquid is intended for use as an assayed quality control material/serum to monitor the precision and trueness of laboratory measurement procedures. All the analyses have been performed by independent collaborative laboratories using SERO's protocol. The mean analytical values and acceptable ranges given in the table below are derived from replicate analyses obtained through collaboration with a number of independent laboratories and they are specific for this LOT.

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The latest data key can be downloaded from our website by using a login.
www.invicon.de

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