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Literatur

- [1] **Aristei C, Kaidar-Person O, Arenas M, Coles C, Offersen BV et al.** (2016): *The 2016 Assisi Think Tank Meeting on breast cancer: white paper.* *Breast Cancer Res. Treat.* 160(2): 211 – 221.
- [2] **Donker M, van Tienhoven G, Straver ME, Meijnen P, van de Velde CJH et al.** (2014): *Radiotherapy or surgery of the axilla after a positive sentinel node in breast cancer (EORTC 10981-22023 AMAROS): a randomised, multicentre, open-label, phase 3 non-inferiority trial.* *Lancet Oncol.* 15(12): 1303 – 1310.
- [3] **Curigliano G, Burstein HJ, Winer EP, Gnant, Dubsy P et al.** (2017): *De-escalating and escalating treatments for early-stage breast cancer: the St. Gallen International Expert Consensus Conference on the Primary Therapy of Early Breast Cancer 2017.* *Ann. Oncol.* 28(8): 1700 – 1712.
- [4] **NCCN guidelines Breast Cancer 2017.** https://www.nccn.org/professionals/physician_gls/pdf/breast.pdf (Accessed 1 September, 2017).
- [5] **Giuliano AE, Ballman K, McCall L, Beitsch P, Whitworth PW et al.** (2016): *Locoregional Recurrence After Sentinel Lymph Node Dissection With or Without Axillary Dissection in Patients With Sentinel Lymph Node Metastases: Long-term Follow-up From the American College of Surgeons Oncology Group (Alliance) ACOSOG Z0011 Randomized Trial.* *Ann. Surg.* 264(3): 413 – 420.
- [6] **Jagsi R, Chadha M, Moni J, Ballman K, Laurie F et al.** (2014): *Radiation field design in the ACOSOG Z0011 (Alliance) Trial.* *J Clin Oncol.* 2014, 32(32): 3600 – 3606.
- [7] **Peg V, Espinosa-Bravo M, Vieites B, Vilardell F, Antúnez JR et al.** (2013): *Intraoperative molecular analysis of total tumor load in sentinel lymph node: a new predictor of axillary status in early breast cancer patients.* *Breast Cancer Res. Treat.* 139(1): 87 – 93.
- [8] **Ohi Y, Umekita Y, Sagara Y, Rai Y, Yotsumoto D et al.** (2012): *Whole sentinel lymph node analysis by a molecular assay predicts axillary node status in breast cancer.* *Br. J. Cancer.* 107(8): 1239 – 43.
- [9] **Peg V, Sansano I, Vieites B, Bernet L, Cano R et al.** (2017): *Role of total tumour load of sentinel lymph node on survival in early breast cancer patients.* *Breast.* 33:8 – 13.
- [10] **Algara López MI.** (2016): *The axilla – less surgery, more radiotherapy? In Kirby A (Chair), Regional nodal irradiation for breast cancer. Symposium conducted at the ESTRO 35 conference, Turin, Italy.*
- [11] **OPTimizing irradiation through molecular assessment of lymph node (OPTIMAL).** <https://clinicaltrials.gov/ct2/show/NCT02335957> (Accessed 1 September, 2017).
- [12] **Algara López MI, Sanz J et al.** (2017): *Irradiazione incidentale versus intenzionale negli stadi iniziali del cancro della mammella e con basso carico tumorale linfonodale. Risultati preliminari dall'OPTIMAL trial. Poster presented at the AIS conference, Florence, Italy.*

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Vertrieb Deutschland: Sysmex Deutschland GmbH
Vertrieb Schweiz: Sysmex Suisse AG
Vertrieb Österreich: Sysmex Austria GmbH
EU Bevollmächtigter: Sysmex Europe GmbH
Hersteller: Sysmex Corporation