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The FUJI DRI-CHEM IMMUNO AG Series is a reliable system with a specificity of 100%, which is crucial for rapid diagnostic testing, for the detection of both Influenza A&B and RSV, the major pathogens responsible for outbreaks. The performance of the device shows higher sensitivity for Influenza A, whatever the subtype, compared to Influenza B and RSV. In this retrospective study, the sensitivity of FUJI IMMUNO AG1 seems to be high enough for rapid diagnostic testing. These preliminary retrospective results need to be completed with prospective data. Conclusion: From this study, FUJI DRI-CHEM IMMUNO AG Influenza AB and RSV kits can be of real value in the management of patients during outbreaks such as hospitals and nursing homes due to the following features : 1/ Satisfying sensitivity results - 2/ Easy and simple use of the system including convenient reading of the results.

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Evaluation of the FUJI DRI-CHEM IMMUNO AG FluAB kit and RSV/ADV kit on Clinical Samples by Use of IMMUNO AG1 Analyser

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INTRODUCTION & PURPOSE

In France, the surveillance of acute respiratory infection in the general population is based on clinical samples performed by general practitioners and paediatricians from the "Réseau Unique". Nasal or nasopharyngeal swabs are collected by the doctors from patients presenting an influenza like illness. Samples are sent to the laboratory by post. The clinical samples arrive at the NIC within an average delay of 4 days. In the laboratory, diagnostic methods include influenza detection by real-time reverse transcription-polymerase chain reaction (RT-PCR) and cell culture. Upon arrival at the NIC, each sample is processed for viral RNA extraction on NucliSens easyMAG instrument (bioMérieux) and tested for the detection of Influenza viruses and RSV by real time RT-PCR and subsequently stored at -80 C.

Although PCR is becoming a method of choice for the diagnostic of respiratory viral infections, easy to handle rapid diagnostic tests are used in laboratories, as well as point of care. Here we evaluated the **FUJI DRI-CHEM IMMUNO AG1** (FUJIFILM, Tokyo, Japan) silver amplification immunochromatography assay, for the detection of Influenza A&B and RSV. This system uses photographic development technology to increase the sensitivity of the conventional immunochromatographic assay.

MATERIAL & METHOD

PROTOCOL


1. This **retrospective study** compared the performance of **FUJI DRI-CHEM IMMUNO AG FluAB kit** and **RSV/Adeno kit** with PCR used at the NIC for virological diagnostic. The Influenza RT-PCR diagnostic was an in-house one step real-time RT-PCR targeting the M gene for Influenza A and NS gene for Influenza B and for RSV the RT-PCR detection was done with a commercialized kit (MWS r-gene™, BioMérieux, Lyon, France) both on the Applied Biosystems 7500 platform.

RT-PCR was performed on fresh specimen and the copy number of viral RNA was available for each positive virus detection.

The Influenza clinical samples were chosen from 2 consecutive winter surveillance 2014/15 for the selection of A H3N2 and B-Yamagata-lineage viruses and 2015/16 for the A(H1N1)pdm09 and B-Victoria lineage specimens. The RSV samples were collected during the 2015/16. The negative came from both winter season. The evaluation included 340 nasopharyngeal specimens including 50 positives for each type/subtype of Influenza A and B viruses, 100 RSV and 40 negatives.

The evaluation of the sensitivity of the **FUJI DRI-CHEM IMMUNO AG1** device was performed on frozen selected clinical samples.

2. A **sensitivity evaluation of FUJI FluAB** was conducted with **4 Influenza strains**. The cell cultures of **A/H3N2/Texas/50/2012**, **A/H1N1/California/7/2009**, **B/Brisbane/60/2008-Victoria** lineage and **B/Massachusetts/2/2012-Yamagata** lineage were adjusted to 10⁸ TCID50/ml and tested twice.

RESULTS

1. Out of the 100 positive specimens per virus the **IMMUNO AG** allowed the detection of **46 Influenza A**, **64 Influenza B** and **67 RSV** and all **negative** clinical samples were negative. We determined the **FUJI IMMUNO AG** the detection limit according to the **Ct value** acquired on RT-PCR for Influenza A per subtype, Influenza B per lineage and for RSV.

The comparison results between RT-PCR and **FUJI DRI-CHEM IMMUNO AG** for **Influenza & RSV** are shown in the table below.

	FUJI A+ / PCR A H3N2+	FUJI A- / PCR A H3N2+	FUJI A+ / PCR A H1N1+	FUJI A- / PCR A H1N1+	FUJI B+ / PCR B Vic+	FUJI B- / PCR B Vic+	FUJI B+ / PCR B Yam+	FUJI B- / PCR B Yam+	FUJI RSV+ / PCR RSV+	FUJI RSV- / PCR RSV+
FUJI / PCR	23	27	23	27	27	23	37	13	67	33
Median Ct	24.61	26.79	25.36	26.9	21.61	22.6	19.66	22.23	19.32	28.31
Extreme	[14.45-28.15]	[23.41-32.08]	[17.2-28.8]	[24.1-28.77]	[19-23.82]	[20.16-24.02]	[14.1-22.62]	[20.3-23.56]	[15.2-27.1]	[18-37.67]

2. The **sensitivity level** of **FUJI FluAB** varies from **10⁴ TCID50/ml** for **Influenza A** to **10⁵ – 10⁶ TCID50/ml** for **Influenza B**. The **FUJI FluAB** sensitivity is comparable to other Rapid Diagnostic Tests with analyser previously evaluated.

Type/Sub-type	A H1N1pdm09	A H3N2	B Victoria	B Yamagata
Strains	A/California/7/2009	A/Texas/50/2012	B/Brisbane/60/2008	B/Massachusetts/2/2012
TCID50/ml	10⁴ / 10⁴	10⁴ / 10⁴	10⁵ / 10⁵	10⁵ / 10⁵

DISCUSSION

The **FUJI DRI-CHEM IMMUNO AG1** is a reliable system with a **specificity of 100%**, which is crucial for rapid diagnostic testing, for the detection of both **Influenza A&B** and **RSV**, the major pathogens responsible for outbreaks. The performance of the device shows higher sensitivity for Influenza A, whatever the subtype, compared to Influenza B and RSV. In this retrospective study, the sensitivity of **FUJI IMMUNO AG1** seems to be high enough for rapid diagnostic testing. These preliminary retrospective results need to be completed with prospective data.

Conclusion: From this study, **FUJI DRI-CHEM IMMUNO AG Influenza AB and RSV kits** can be of real value in the management of patients during outbreaks such as hospitals and nursing homes due to the following features : 1/ Satisfying sensitivity results - 2/ Easy and simple use of the system including convenient reading of the results.

REFERENCES

- M. Bouscambert Duchamp, J. S. Casalegno, Y. Gillet, E. Frobert, E. Bernard, V. Escuret, G. Billaud, M. Valette, E. Javouhey, B. Lina, D. Floret, F. Morfin. Pandemic A (H1N1) 2009 influenza virus detection by real time RT-PCR: is viral quantification useful? Clin Microbiol Infect. 2010; 16: 317-321
- K. Mitamura, H. Shimizu, M. Yamazaki, M. Ichikawa, K. Nagai, J. Katada, A. Wada, C. Kawakami, N. Sugaya. Clinical evaluation of highly sensitive silver amplification immunochromatography systems for rapid diagnosis of influenza. J. Virol. Method. 2013;194:123-128

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Niederlassung der FUJIFILM Europe GmbH