

# COVID-19 Neutralizing Antibody (NAb)


WHO-Standard:  
Einheit in BAU/ml



Based on Neutralizing Antibody (NAb) level, the Titer is displayed as a Numeric Value in BAU/ml followed by Result with an audible signal:

NAb Titer Range	Result
• < 30 BAU/ml	NEGATIVE
• ≥ 30 - < 250 BAU/ml	LOW +
• ≥ 250 - < 500 BAU/ml	MEDIUM +
• ≥ 500 BAU/ml	HIGH +

## Testdurchführung Antikörper Bestimmung:



Bei Interesse wenden Sie sich gerne an folgenden Kontakt:

**Telefon: +43 1 486 16 31**

**E-Mail: [aktionen@sysmex.at](mailto:aktionen@sysmex.at)**